



## APPLICATION NOTE

Quick DNA Extraction from Rice Seed (Wet)

With kind permission of RiceTec Inc, Alvin TX

## ABSTRACT

DNA extractions can be a very time consuming and tedious process. Finding a quick method in which DNA could be extracted and used for PCR is essential. Described below is a quick "dirty" method that produces a high enough concentration of DNA that can be used for PCR.

## Sample Extraction

Samples are prepared using a 96 well 1ml assay block. Dispense one 5/32" (4mm) stainless steel bead into each well using the Grinding ball dispenser (SPEX SamplePrep cat. # 2100). Next, add one seed to each well. Dispense extraction buffer into each well and securely cap each well. After the samples have been capped, grind them in the Geno/Grinder at 500 strokes/minute for two minutes. Centrifuge for 1 min to bring all liquid to the bottom of the assay block. Incubate the samples in about 1inch (25 mm) of water at 95°C for 20 minutes then place them on ice for approximately 10 minutes or until samples are cool to the touch. Centrifuge again for 1 minute. Add neutralizing extraction buffer and seal the assay block with sealing film. Centrifuge the samples for 10 minutes at 3000rpm. Transfer 300µL of the supernatant to a clean 96 well plate. DNA can be further purified with clean-up kits available on the market.

## Yields

The total concentrations yielded from the samples range from 3-7ng/µL in a final volume of 200µL with purities of 1.5- 1.8 ng/µL.

## Conclusion

The method described above is sufficient for PCR and it takes less time than the standard chloroform extraction. Total time ranges from one to two hours.

**:: APPLICATION NOTE SP025:**  
Cell Disruption

**:: APPARATUS:**  
**Geno/Grinder®**

**:: APPLICATION:**  
DNA Extraction from Rice  
Seeds



### SPEX SamplePrep

65 Liberty St  
Metuchen, NJ 08840 USA  
Tel: 732-623-0465  
Fax: 732-906-2492

E-mail: [Sampleprep@spex.com](mailto:Sampleprep@spex.com)  
[www.spexsampleprep.com](http://www.spexsampleprep.com)

### European Headquarters

SPEX Europe  
2 Dalston Gardens  
Stanmore, HA7 1BQ, UK  
Tel: +44 (0)208 204 6656  
Fax: +44 (0)208 204 6654  
E-mail: [spexeurope@spex.com](mailto:spexeurope@spex.com)  
Web: [www.spexeurope.com](http://www.spexeurope.com)