


Introduction

Decommissioning of nuclear reactor sites requires accurate characterisation of multiple radionuclides in a range of complex sample matrices. Quantification of difficult-to-measure alpha- and beta-emitters requires complete digestion of samples prior to radiochemical separation and measurement. The Katanax K2 has been applied to a number of decommissioning matrices, allowing automated digestion of samples by lithium borate and alkali fusion without the need for handling of hazardous reagents e.g. hydrofluoric acid.



Applications


Concrete

- Activation product ^{41}Ca
- 0.5 g blank concrete
- Flux composition 1:1:0.25 lithium metaborate : lithium tetraborate : lithium bromide wetting agent
- Melt dissolved in 5M HNO_3 followed by PEG precipitation
- 97% Ca recovered after the PEG precipitation.



Graphite

- Multiple radionuclides



Sand, tuff and titanium dioxide

- Reference materials characterised for naturally occurring radioactivity content (U and Th decay series)
- Fusion flux varies depending on sample matrix:
 - **Sand and tuff:** 1:1 lithium metaborate : lithium tetraborate
 - **TiO₂:** alkali fusion with 4:3 Na_2CO_3 : H_3BO_3

